

Rec'd PCT/PTO 03 JUN 2005



10/557435

GB 03 / 5272

PRIORITY DOCUMENT
SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH
RULE 17.1(a) OR (b)

The Patent Office
Concept House
Cardiff Road
Newport
South Wales
NP10 8QQ

REC'D 02 MAR 2004

WIPO

PCT

I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation & Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.

In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.

In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.

Signed

Stephen Hordley

Dated

22 January 2004


Best Available Copy



1/77

Request for grant of a patent

(See the notes on the back of this form. You can also get an explanatory leaflet from the Patent Office to help you fill in this form)



03DEC02 E768074-5 D02884
P01/7700 0.00-0228139.2

The Patent Office

Cardiff Road
Newport
South Wales
NP10 8QQ

1. Your reference

P32365/LMC/GST

2. Patent application number

(The Patent Office will fill in this part)

0228139.2

03 DEC 2002

3. Full name, address and postcode of the or of each applicant (underline all surnames)

Dr Olga Kozlova-Zwiderman
239/5 Gilmerton Road
Edinburgh
EH16 5TH

Patents ADP number (if you know it)

If the applicant is a corporate body, give the country/state of its incorporation

08519084001

4. Title of the invention

"Fungal Biosensor for Contaminant Detection"

5. Name of your agent (if you have one)

Murgitroyd & Company

"Address for service" in the United Kingdom to which all correspondence should be sent (including the postcode)

Scotland House
165-169 Scotland Street
Glasgow
G5 8PL

Patents ADP number (if you know it)

1198013

6. If you are declaring priority from one or more earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (if you know it) the or each application number

Country

Priority application number
(if you know it)

Date of filing
(day / month / year)

7. If this application is divided or otherwise derived from an earlier UK application, give the number and the filing date of the earlier application

Number of earlier application

Date of filing
(day / month / year)

8. Is a statement of inventorship and of right to grant of a patent required in support of this request? (Answer 'Yes' if:

No

- a) any applicant named in part 3 is not an inventor, or
 - b) there is an inventor who is not named as an applicant, or
 - c) any named applicant is a corporate body.
- See note (d))

Patents Form 1/77

Enter the number of sheets for any of the following items you are filing with this form. Do not count copies of the same document

Continuation sheets of this form

Description

19

Claim(s)

-

Abstract

-

Drawing(s)

9

19 ph

10. If you are also filing any of the following, state how many against each item.

Priority documents

-

Translations of priority documents

-

Statement of inventorship and right to grant of a patent (Patents Form 7/77)

-

Request for preliminary examination and search (Patents Form 9/77)

-

Request for substantive examination (Patents Form 10/77)

-

Any other documents (please specify)

-

11. I/We request the grant of a patent on the basis of this application.

Signature Murgitroyd & Company

Date

Murgitroyd & Company

2 December 2002

12. Name and daytime telephone number of person to contact in the United Kingdom

Gordon Stark

0141 307 8400

Warning

After an application for a patent has been filed, the Comptroller of the Patent Office will consider whether publication or communication of the invention should be prohibited or restricted under Section 22 of the Patents Act 1977. You will be informed if it is necessary to prohibit or restrict your invention in this way. Furthermore, if you live in the United Kingdom, Section 23 of the Patents Act 1977 stops you from applying for a patent abroad without first getting written permission from the Patent Office unless an application has been filed at least 6 weeks beforehand in the United Kingdom for a patent for the same invention and either no direction prohibiting publication or communication has been given, or any such direction has been revoked.

Notes

- If you need help to fill in this form or you have any questions, please contact the Patent Office on 08459 500505.
- Write your answers in capital letters using black ink or you may type them.
- If there is not enough space for all the relevant details on any part of this form, please continue on a separate sheet of paper and write "see continuation sheet" in the relevant part(s). Any continuation sheet should be attached to this form.
- If you have answered 'Yes' Patents Form 7/77 will need to be filed.
- Once you have filled in the form you must remember to sign and date it.
- For details of the fee and ways to pay please contact the Patent Office.

1 **"Fungal Biosensor"**

2

3 The present invention provides a method of using
4 fungus for determining the presence of at least one
5 toxic substance in a sample and for assisting in the
6 identification of the toxicant(s). More
7 specifically there is provided a toxicity assay for
8 use in determining the presence of toxins and in
9 particular heavy metals and organophenols.

10

11 The release of contaminating substances into an
12 environment such as a waterway or an area of
13 agricultural land can have serious effects on the
14 ecosystems found in that environment. It is
15 important to be able to analyse these effects both
16 prior to the release of such contaminants so as to
17 manage their treatment or release, and after release
18 so as to determine and counteract their effects.

19

20 Current methods used to monitor water quality and
21 screen effluent generally involve chemical toxicity
22 tests. However, these tests require a general idea

1 of the type of contaminant being tested for and can
2 be very expensive.

3

4 Biosensors are also used for toxicity testing and
5 are well known in the field. Toxicity depends on a
6 variety of factors including pH, temperature,
7 salinity and contaminant concentration, but depends
8 especially on the test organism used in the sensor.

9

10 One of the most commonly used organisms is the
11 bioluminescent bacterium, *Vibrio fischeri*. The
12 bioluminescence involved is mediated by the
13 luciferin-luciferase enzyme system wherein light
14 emission is dependent on the electron transfer
15 chain. Any disruption to the electron transfer
16 chain, for example on exposure to a toxicant,
17 affects light emission. Measurement of changes in
18 light emission is therefore indicative of the
19 presence of a toxic substance.

20

21 This system, however, only provides a simple
22 indication of whether a contaminant is toxic or not.
23 No detailed information is obtained on how toxic the
24 contaminant is, nor is the contaminant identified.

25

26 The term toxicant as herein described relates to
27 compounds which are toxic to fungus and also to any
28 substance or compound with anti-fungal activity.

29

30 According to the present invention there is provided

- 1 an assay for use in determining the presence of a
2 known toxicant in a test sample, the assay
3 comprising the steps of;
- 4 - exposing a fungi transformed with a
5 recombinant aequorin gene to a test sample of
6 a substance,
 - 7
 - 8 - measuring the luminescence produced by the
9 fungi,
 - 10
 - 11 - converting the luminescence data into a
12 cytosolic free calcium ion concentration
13 trace,
 - 14
 - 15 - and comparing at least one parameter of the
16 cytosolic free calcium ion concentration
17 trace with a bank of known toxicity reference
18 data.

19
20 Preferably the cytosolic free calcium ion trace is a
21 plot of the cytosolic free calcium ion concentration
22 against time.

23
24 Preferably the fungi transformed with a recombinant
25 aequorin gene is a filamentous fungi.

26
27 More preferably the fungi is of the *Aspergillus*
28 species.

29
30 Most preferably the fungi is *Aspergillus awamori*.

31

1 Most preferably still the strain of *Aspergillus*
2 *awamori* is strain 66A.

3

4 Preferably the substance is a contaminant.

5

6 Preferably the substance is a contaminated sample.

7

8 Preferably the parameter is at least one or more
9 selected from the group comprising; lag time, rise
10 time, amplitude, full width half maximum, number of
11 cytosolic free calcium ion concentration increases,
12 percentage increase in final cytosolic free calcium
13 ion concentration resting level, and percentage
14 increase in recovery time.

15

16 In a further embodiment of the invention the test
17 sample is added in advance of the application of a
18 stimulus to the test sample.

19

20 Preferably the stimulus is at least one or more from
21 the group comprising; mechanical perturbation, hypo-
22 osmotic shock, and change in external calcium
23 chloride concentration.

24

25 Preferably the test sample is added 1 minute to 1
26 hour prior to the application of the stimulus.

27

28 More preferably the test sample is added 5 minutes
29 prior to the application of the stimulus.

30

31 More preferably the test sample is added 30 minutes
32 prior to the application of the stimulus.

1
2 In such instances, the parameters may include at
3 least one or more selected from the group
4 comprising; lag time, rise time, amplitude, full
5 width half maximum, number of cytosolic free calcium
6 ion concentration increase, percentage increase in
7 final cytosolic free calcium ion concentration
8 resting level, percentage increase in recovery time,
9 and percentage increase in pre-stimulating cytosolic
10 free calcium ion concentration resting level.

11
12 Preferably luminescence is measured for between 1
13 minute and 5 hours following the application of the
14 stimulus.

15
16 More preferably luminescence is measured for 5
17 minutes following the application of the stimulus.

18
19 A further embodiment of the present invention
20 provides a method of determining the amount of a
21 contaminant present in a test sample comprising the
22 steps of;

- 23
- 24 - exposing a strain of fungi transformed with a
 - 25 recombinant aequorin gene to a sample,
 - 26
 - 27 - measuring the luminescence produced by the
 - 28 fungi,
 - 29
 - 30 - converting the luminescence data into a
 - 31 cytosolic free calcium ion concentration
 - 32 trace,

1
2 - and comparing at least one parameter of the
3 cytosolic free calcium ion concentration
4 trace with a bank of pre-prepared toxicity
5 reference data.

6
7 Preferably the cytosolic free calcium ion trace is a
8 plot of the cytosolic free calcium ion concentration
9 against time.

10
11 Preferably the fungi transformed with a recombinant
12 aequorin gene is a filamentous fungi.

13
14 More preferably the fungi is of the *Aspergillus*
15 species.

16
17 Most preferably the fungi is *Aspergillus awamori*.

18
19 Most preferably still the strain of *Aspergillus*
20 *awamori* is strain 66A.

21
22 Preferably the substance is a contaminant.

23
24 Preferably the substance is a contaminated sample.

25
26 Preferably the parameter is at least one or more
27 selected from the group comprising; lag time, rise
28 time, amplitude, full width half maximum, number of
29 cytosolic free calcium ion concentration increases,
30 percentage increase in final cytosolic free calcium
31 ion concentration resting level, and percentage
32 increase in recovery time.

1

2 In a further embodiment of the invention the test
3 sample is added in advance of the application of a
4 stimulus.

5

6 Preferably the stimulus is at least one or more from
7 the group comprising; mechanical perturbation, hypo-
8 osmotic shock, and change in external calcium
9 chloride concentration.

10

11 Preferably the test sample is added 1 minute to 1
12 hour prior to the application of the stimulus.

13

14 More preferably the test sample is added 5 minutes
15 prior to the application of the stimulus.

16

17 More preferably the test sample is added 30 minutes
18 prior to the application of the stimulus.

19

20 In such instances, the parameters may include at
21 least one or more selected from the group
22 comprising; lag time, rise time, amplitude, full
23 width half maximum, number of cytosolic free calcium
24 ion concentration increase, percentage increase in
25 final cytosolic free calcium ion concentration
26 resting level, percentage increase in recovery time,
27 and percentage increase in pre-stimulating cytosolic
28 free calcium ion concentration, resting level.

29

30 Preferably luminescence is measured for between 1
31 minute and 5 hours following the application of the
32 stimulus.

1
2 More preferably luminescence is measured for 5
3 minutes following the application of the stimulus.
4

5 A first experiment comprises testing the effect of
6 pre-incubation of *Aspergillus awamori* with toxicants
7 on cytosolic free calcium ion concentration response
8 to an increase in external calcium chloride.
9

10 A further set of experiments described herein shows
11 attempts to obtain characteristic data for a range
12 of different toxicants at a number of different
13 concentrations. The results demonstrate that each
14 toxicant at each concentration produces a
15 distinctive cytosolic free calcium ion concentration
16 trace whose traits could be used to identify and
17 characterise a toxicant present in a test sample.
18

19 A final experiment attempts to determine whether it
20 is possible to identify and characterise individual
21 toxicants from testing samples of mixtures of
22 toxicants in different proportions. The traces
23 produced are distinct for each mixture.
24

25 These results show that it is possible to
26 characterise and identify a specific toxicant from a
27 test sample by using the characteristic data
28 obtained from a cytosolic free calcium ion
29 concentration trace.
30

31 The parameters referred to herein relate to the
32 following;

- 1
- 2 Lag Time, the time from addition of the test sample
- 3 to the time when the cytosolic free calcium ion
- 4 concentration, $[Ca^{2+}]_c$, began to rise;
- 5
- 6 Rise Time, the time from addition of the test sample
- 7 to the time at which maximum $[Ca^{2+}]_c$ was reached;
- 8
- 9 Amplitude, the maximum $[Ca^{2+}]_c$ reached during the
- 10 experiment;
- 11
- 12 Full Width Half Maximum, the width of the transient
- 13 at the point where the amplitude equals half of the
- 14 maximum amplitude of the transient;
- 15
- 16 Number of $[Ca^{2+}]_c$ Rises, the number of peaks in
- 17 $[Ca^{2+}]_c$;
- 18
- 19 Percentage Increase in Final $[Ca^{2+}]_c$ Resting Level,
- 20 the percentage increase in resting $[Ca^{2+}]_c$ at the end
- 21 of the experiment, where the control value is taken
- 22 to be 100%;
- 23
- 24 Percentage Increase in Recovery Time, percentage
- 25 increase in recovery time where recovery time
- 26 represents the total amount of $[Ca^{2+}]_c$ released
- 27 during the period of time from the point when the
- 28 maximum amplitude following calcium chloride
- 29 treatment was achieved to the point when the $[Ca^{2+}]_c$
- 30 reached its final resting level. Recovery time was
- 31 initially calculated for control cultures. In the
- 32 control this period of time was calculated as 250

seconds. For the cultures subjected to the treatment with toxicant(s) the total amount of $[Ca^{2+}]_c$ was calculated for the same period of 250 seconds starting from the maximum amplitude. The recovery time of the control cultures was therefore:

$$\frac{\text{total amount of } [Ca^{2+}]_c \text{ (}\mu\text{M) for the toxicant-treated samples over 250 seconds} \times 100}{\text{total amount of } [Ca^{2+}]_c \text{ (}\mu\text{M) for the control sample over 250 seconds}}$$

; and

Percentage Increase in pre-Stimulating $[Ca^{2+}]_c$ Resting Level, the percentage increase in $[Ca^{2+}]_c$ prior to the stimulus, where the control value is taken to be 100%.

There is also the possibility of combining one or more of these parameters to obtain further values which can be used for identification of the toxicants in the mixture. For example, the summation of amplitude and recovery time will give the value of total cytosolic free calcium ions emitted from the time when $[Ca^{2+}]_c$ reaches its peak. Also summation of lag time and rise time will give the total time required for $[Ca^{2+}]_c$ to reach its peak. The division of final $[Ca^{2+}]_c$ resting level onto the pre-stimulation $[Ca^{2+}]_c$ resting level will show how many times the $[Ca^{2+}]_c$ resting level has changed after stimulation. Similarly, a division of the final $[Ca^{2+}]_c$ resting level onto the initial

1 $[Ca^{2+}]_c$ resting level prior to the addition of
2 toxicant(s) gives further identifying data.
3 Additionally, the summation of all the data points
4 of the trace gives the total amount of cytosolic
5 free calcium ions released during the monitoring
6 period.

7
8 The present invention will now be described with
9 reference to the following non-limiting examples and
10 with reference to the figures, wherein:

11
12 Figure 1 shows the characteristic $[Ca^{2+}]_c$ trace
13 produced on addition of 5mM external $CaCl_2$,
14 following a 5 minute pre-incubation with
15 different concentrations of 3,5-DCP.

16
17 Figure 2 shows the characteristic $[Ca^{2+}]_c$ trace
18 produced on addition of 5mM external $CaCl_2$,
19 following a 5 minute pre-incubation with
20 different concentrations of Cr^{6+} .

21
22 Figure 3 shows the characteristic $[Ca^{2+}]_c$ trace
23 produced on addition of 5mM external $CaCl_2$,
24 following a 5 minute pre-incubation with
25 different concentrations of Zn^{2+} .

26
27 Figure 4 shows the characteristic $[Ca^{2+}]_c$ trace
28 produced on addition of 5mM external $CaCl_2$,
29 following a 30 minute pre-incubation with
30 different concentrations of 3,5-DCP.

31

1 Figure 5 shows the characteristic $[Ca^{2+}]_c$ trace
2 produced on addition of 5mM external $CaCl_2$,
3 following a 30 minute pre-incubation with
4 different concentrations of Cr^{6+} .

5
6 Figure 6 shows the characteristic $[Ca^{2+}]_c$ trace
7 produced on addition of 5mM external $CaCl_2$,
8 following a 30 minute pre-incubation with
9 different concentrations of Zn^{2+} .

10
11 Figure 7 shows the characteristic cytosolic free
12 calcium ion concentration, $[Ca^{2+}]_c$, trace
13 produced on addition of 5mM $CaCl_2$ following a 5
14 minute pre-incubation with different
15 concentrations of 3,5-dichlorophenol, 3,5-DCP.

16
17 Figure 8 shows the characteristic $[Ca^{2+}]_c$ trace
18 produced on addition of 5mM $CaCl_2$, following a
19 30 minute pre-incubation with different
20 concentrations of 3,5-DCP.

21
22 Figure 9 shows the characteristic $[Ca^{2+}]_c$ trace
23 produced on addition of 5mM $CaCl_2$, following a 5
24 minute pre-incubation with different
25 concentrations of chromium ions, Cr^{6+} .

26
27 Figure 10 shows the characteristic $[Ca^{2+}]_c$ trace
28 produced on addition of 5mM $CaCl_2$, following a
29 30 minute pre-incubation with different
30 concentrations of chromium ions, Cr^{6+} .

31

1 Figure 11 shows the characteristic $[Ca^{2+}]_c$ trace
2 produced on addition of 5mM $CaCl_2$, following a 5
3 minute pre-incubation with different
4 concentrations of zinc ions, Zn^{2+} .

5
6 Figure 12 shows the characteristic $[Ca^{2+}]_c$ trace
7 produced on addition of 5mM $CaCl_2$, following a
8 30 minute pre-incubation with different
9 concentrations of zinc ions, Zn^{2+} .

10
11 Figure 13 shows the values of $[Ca^{2+}]_c$ trace
12 parameters characteristic for different
13 concentrations of pentochlorophenol, PCP; sodium
14 dodecyl sulphate, SDS; and Toluene. Parameters
15 assessed are Lag Time, LT; Rise Time, RT;
16 Amplitude, A; Full Width Half Maximum, FWHM;
17 Percentage Increase in pre-Stimulating $[Ca^{2+}]_c$.
18 Resting Level, %IpreSRL; Percentage Increase in
19 Final $[Ca^{2+}]_c$ Resting Level, %IFRL; Percentage
20 Increase in Recovery Time, %IRT; and Number of
21 $[Ca^{2+}]_c$ Increases.

22
23 Figure 14 shows the values of $[Ca^{2+}]_c$ trace
24 parameters characteristic for 3,5-DCP, PCP, Zn^{2+} ,
25 Cr^{6+} , Toluene, and SDS. Parameters assessed are
26 Lag Time, LT; Rise Time, RT; Amplitude, A; Full
27 Width Half Maximum, FWHM; Percentage Increase in
28 pre-Stimulating $[Ca^{2+}]_c$ Resting Level, %IpreSRL;
29 Percentage Increase in Final $[Ca^{2+}]_c$ Resting
30 Level, %IFRL; Percentage Increase in Recovery
31 Time, %IRT; and Number of $[Ca^{2+}]_c$ Increases.

32

Figure 15 shows the values of $[Ca^{2+}]_e$ trace parameters characteristic for different mixtures of toxicants. Parameters assessed are Lag Time, LT; Rise Time, RT; Amplitude, A; Full Width Half Maximum, FWHM; Percentage Increase in pre-Stimulating $[Ca^{2+}]_e$ Resting Level, %IpreSRL; Percentage Increase in Final $[Ca^{2+}]_e$ Resting Level, %IFRL; Percentage Increase in Recovery Time, %IRT; and Number of $[Ca^{2+}]_e$ Increases.

Effect of pre-incubation of *Aspergillus awamori* with toxicants on $[Ca^{2+}]_e$ response to external calcium chloride

12 ml of sterile VS medium was inoculated with 1×10^5 spores per ml *A. awamori* strain 66A. 100 μ l of the inoculated medium was added to each well of a 96-well plate and cultured in a humidity chamber in the presence of free water at 30 °C for 24 hours.

The following toxicants were tested: 3,5-dichlorophenol, zinc sulphate, and potassium dichromate. Each toxicant was added in a total volume of 25 μ l VS medium or water 5 or 30 minutes before addition of 5 mM calcium chloride.

Luminescence was monitored for 5 minutes following addition of $CaCl_2$. Aequorin was completely discharged by adding 3M calcium chloride in 20% ethanol. The total concentration is thus 1.5 M calcium chloride in 10% ethanol.

1 Luminometry was performed using an EG & G Berthold
2 (Bad Wildbad, Germany) LB96P Microlumat luminometer.
3 Luminescence data was converted from real light
4 units to $[Ca^{2+}]_i$ values using the following equation:

5
6
$$PCa = 0.332588 (-\log k) + 5.5593,$$

7
8 where k = luminescence counts per second/total
9 luminescence counts. Total luminescence is measured
10 as an integral of all luminescence up to complete
11 aequorin discharge.

12

13 The following parameters were assessed:

14 Rise Time, Amplitude, Full Width Half Maximum and
15 Final $[Ca^{2+}]_i$ Resting Level.

16

17 Effects of different concentrations of toxicants on
18 $[Ca^{2+}]_i$ traces

19

20 *Aspergillus awamori* were transformed with an
21 expression vector (pAEQ1-15) comprising a gene for
22 synthetic apoaequorin (*aeqS*) under the control of
23 the constitutive glucose-6-phosphate dehydrogenase
24 promoter (*gpdA*).

25

26 These transformants were cultured in 100 μ l of
27 Vogel's medium with 1% sucrose (VS medium) in
28 microwell plates for 24 hours before addition of a
29 toxicant or a control of distilled water. Toxicants
30 were dissolved in water to give the concentrations
31 shown below. 25 μ l of the each of the following
32 concentrations were added to each culture:

1

TOXICANT	CONCENTRATIONS (mg/l)
3,5-dichlorophenol (3,5-DCP)	0.112, 11.2, 112
Chromium ions (Cr^{6+})	15, 120, 260
Zinc ions (Zn^{2+})	180, 350, 700, 1300

2

3 The cultures were incubated for 5 or 30 minutes
 4 before addition of 100 μl 5mM CaCl_2 . Luminescence
 5 was measured for 5 minutes using a plate
 6 luminometer. Luminescence data was manually
 7 converted from relative light units to cytosolic
 8 free calcium ion concentration, $[\text{Ca}^{2+}]_c$. This was
 9 then plotted against time and parameters of this
 10 trace were analysed. Parameters assessed were as
 11 follows:

12

13 Rise Time, the time from addition of CaCl_2 to the
 14 moment when maximum $[\text{Ca}^{2+}]_c$ was achieved;

15 Amplitude, the maximum $[\text{Ca}^{2+}]_c$ reached during the
 16 experiment;

17 Full Width Half Maximum, the width of the transient
 18 at the point where the amplitude equals half of the
 19 maximum amplitude of the transient;

20 and Final Resting $[\text{Ca}^{2+}]_c$ Level, the resting $[\text{Ca}^{2+}]_c$
 21 at the end of the experiment.

22

23 Effects of further toxicants on $[\text{Ca}^{2+}]_c$ traces

24

25 Cultures of *Aspergillus awamori* as described above
 26 were used to test the effects of further toxicants.
 27 The concentrations of toxicants tested were made up

1 as follows in water, where the concentrations tested
2 were based on Dutch target and intervention values
3 for toxicants and Kelly Guidelines for the
4 classification of contaminated soils:
5

TOXICANT	CONCENTRATION (mg/l)
Pentachlorophenol, PCP	0.01, 0.1, 1, 5, 10
Sodium dodecyl sulphate, SDS	1, 10, 50, 100, 500
Toluene	1, 25
3,5-DCP	10
Zn ²⁺	700
Cr ⁶⁺	15

6
7
8 In the first set-up (S1), 100 µl of each toxicant
9 concentration or of the control (VS medium) were
10 added to the cultures through built-in injectors and
11 luminescence monitored for 5 minutes. In a second
12 set of experiments (S2), cultures were pre-incubated
13 with the toxicant or control for 5 minutes before
14 addition of 5mM CaCl₂ in a total volume of 25 µl
15 distilled water. Luminescence was monitored for 5
16 minutes following addition of CaCl₂. Luminescence
17 data was converted from relative light units to
18 [Ca²⁺]_c values as described above. The following
19 parameters were assessed in S1:
20 Lag Time, the time from addition of CaCl₂ to the
21 time when [Ca²⁺]_c began to rise;
22 Rise Time;
23 Amplitude;
24 Full Width Half Maximum;

1 Percentage Increase in Final $[Ca^{2+}]_0$ Resting Level,
2 where the control value was taken to be 100%;
3 Percentage Increase in Recovery Time, where the
4 control value was taken to be 100%;
5 and Number of $[Ca^{2+}]_0$ Increases, the number of $[Ca^{2+}]_0$
6 transients.

7
8 In S2, the Percentage Increase in pre-Stimulating
9 $[Ca^{2+}]_0$ Resting Level, where the control value was
10 taken to be 100%, was assessed in addition to all of
11 the parameters tested in S1.

12
13 Effects of mixtures containing different proportions
14 of toxicants on $[Ca^{2+}]_0$ traces

15
16 The experiments described when examining the effects
17 of further toxicants were repeated for different
18 mixtures of toxicants. The following mixtures were
19 made up in water for testing:

20
21 6 mg/l 3,5-DCP + 12 mg/l Cr^{6+}
22 30 mg/l Cr^{6+} + 350 mg/l Zn^{2+}
23 10 mg/l 3,5-DCP + 350 mg/l Zn^{2+}
24 6 mg/l 3,5-DCP + 12 mg/l Cr^{6+} + 350 mg/l Zn^{2+}

25 Mixture 1: 20 mg/l Cadmium
26 100 mg/l Copper
27 50 mg/l Chromium
28 250 mg/l Zinc
29 500 mg/l SDS

30 Mixture 2: 20 mg/l Cadmium
31 100 mg/l Copper
32 50 mg/l Chromium

250 mg/l Zinc

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

These experiments demonstrate a novel finding that each toxicant results in a different and characteristic $[Ca^{2+}]_c$ transient. Additionally each concentration of toxicant produces a unique $[Ca^{2+}]_c$ transient. From these characteristic fingerprint responses a profile of data can be built up and used to create a bank of data for each toxicant. Results from testing samples can be compared with this data bank and the presence of a particular toxicant can thus be determined. Furthermore, details such as the mode of action of the toxicant, and the amount of toxicant present can be deduced from a comparison with the bank of pre-gathered data.

Although the invention has been particularly shown and described with reference to particular examples, it will be understood by those skilled in the art that various changes in the form and details may be made therein without departing from the scope of the present invention.




Figure 1




Figure 2




Figure 3




Figure 4




Figure 5




Figure 6




Figure 7




Figure 8




Figure 9




Figure 10




Figure 11




Figure 12

Chemical	Values of interest (mg/l)	LT		RT		A		FWHM		%IpreSRL		%IFRL		%IRT		Number of increases	
		S1	S2	S1	S2	S1	S2	S1	S2	S2	S2	S1	S2	S1	S2	S1	S2
PCP	0.01	0	0	1	1	75±23	100±5	↑	-	97±9	104±6	115±6	104±6	111±4	102±3	1	1
	0.1	0	0	1	1	74±22	91±5	↑	-	104±9	104±11	118±5	104±11	122±9	99±3	1	1
	1	0	0	1	1	72±14	80±10	↑	-	109±5	107±11	204±13	107±11	195±8	109±6	1	1
	5	0	0	12.6	1	199±26	89±14	↑	↑	120±3	205±9	209±7	205±9	274±4	218±7	1	1
	10	0	0	12.6	12.6	417±73	175±27	↑	↑	131±9	253±7	305±17	253±7	373±22	308±4	1	1
SDS	1	0	0	1	1	111±48	102±8	-	-	125±34	106±7	116±5	106±7	114±6	109±9	1	1
	10	0	0	12.6	1	246±19	118±4	↑	↑	136±13	120±13	162±7	120±13	243±14	154±10	1	1
	50	0	0	12.6	1	295±33	115±12	↑	↑	323±19	222±16	237±9	222±16	328±6	287±10	1	1
	100	0	0	12.6	1	293±19	116±3	↑	↑	359±24	256±14	405±15	256±14	367±11	286±12	1	1
	500	0	n.a.	12.6	n.a.	998±35	n.a.	↑	n.a.	n.a.	n.a.	561±12	n.a.	565±12	n.a.	1	n.a.
Toluene	1	n.a.	0	n.a.	1	n.a.	74±13	n.a.	-	129±9	113±3	n.a.	113±3	n.a.	115±2	n.a.	1

Note: LT= lag time

RT=rise time

A=changes in A (%)

FWHM=full width half maximum

%IpreSRL= % increase in pre-stimulation resting level

%IFRL=% increase in final resting level

%IRT=% increase in recovery time

S1=Stage 1

S2=Stage 2

Figure 13

Chemical	Values of interest (mg/l)	LT		RT		A		FWHM	%IpreSRL	%IFRL		%IRT		Number of increases		
		S1	S2	S1	S2	S1	S2			S1	S2	S1	S2			
3,5 DCP	10	0	0	12.6	1	40±2	79±25	↑	↑	102±3	134±10	172±21	153±13	207±23	1	1
PCP	10	0	0	12.6	12.6	417±73	175±27	↑	↑	131±9	305±17	253±7	373±22	308±4	1	1
Zn ²⁺	700	0	0	12.6	1	42±5	74±1	↑	↑	142±18	225±8	263±14	258±5	225±12	1	1
Cr ⁶⁺	15	0	0	1	1	41±5	84±23	-	-	102±12	104±5	110±11	103±2	120±28	1	1
Toluene	25 (1)	0	0	n.a.	1	n.a.	73±13	n.a.	↑	121±6	n.a.	113±3	n.a.	115±2	1	1
SDS	500	0	0	12.6	1	998±35	116±3	↑	↑	359±24	561±12	256±14	565±12	286±12	1	1

Note: LT= lag time

RT=rise time

A=changes in A (%)

FWHM=full width half maximum

%IpreSRL= % increase in pre-stimulation resting level

%IFRL= % increase in final resting level

%IRT= % increase in recovery time

S1=Stage 1

S2=Stage 2

Figure 14

Chemical	Values of interest (mg/l)		LT		RT		A		FWHM		%IpreSRL		%IFRL		%IRT		Number of increases	
	S1	S2	S1	S2	S1	S2	S1	S2	S1	S2	S1	S2	S1	S2	S1	S2	S1	S2
3,5 DCP	0	0	12.6	1	40±2	79±25	↑	↑	102±3	134±10	172±21	153±13	207±23	1	1	1	1	
Cr ⁶⁺	0	0	1	1	41±5	84±23	-	-	102±12	104±5	110±11	103±2	120±28	1	1	1	1	
Zn ²⁺	0	0	12.6	1	42±5	74±1	↑	↑	142±18	225±8	263±14	258±5	225±12	1	1	1	1	
SDS	0	0	12.6	1	998±35	116±3	↑	↑	359±24	561±12	256±14	565±12	286±12	1	1	1	1	
3,5-DCP + Cr ⁶⁺	0	0	12.6	1	33±4	88±13	↑	-	96±26	117±4	120±14	128±5	119±10	1	1	1	1	
Cr ⁶⁺ Zn ²⁺	0	0	1	1	23±6	76±4	↑	-	100±5	143±18	150±12	158±26	154±6	1	2	1	2	
3,5DCP + Zn ²⁺	0	0	12.6	1	65±5	86±2	↑	-	103±9	294±18	153±16	311±14	208±11	1	2	1	2	
3,5-DCP + Cr ⁶⁺ Zn ²⁺	0	0	12.6	1	25±2	79±5	↑	-	102±8	164±4	150±9	160±5	195±8	1	2	1	2	
Mixture 1	0	0	12.6	12.6	466±13	128±6	↑	↑	262±13	402±17	501±38	446±17	477±28	2	2	2	2	
Mixture 2	0	0	1	1	116±8	69±14	↑	-	148±16	177±44	132±27	170±28	120±7	1	1	1	1	

Italics represents data obtained with very high concentrations of toxicants: Zn²⁺ = 700 mg/l; Cr⁶⁺ = 120 mg/l; 3,5 DCP = 49 mg/l

Figure 15

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☒ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☒ GRAY SCALE DOCUMENTS
- ☒ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.